


InLab[®] Redox Micro

Micro ORP electrode with platinum ring

Scale 1:1



Sensor type combined ORP electrode

Order number 51343203

Specifications

pH range -
zero point (pH) -

temperature range 0...100 °C
temperature probe -

shaft material glass
membrane glass -
membrane resistance -

reference system ARGENTHAL™
reference diaphragm ceramic
reference electrolyte 3 mol/L KCl
bridge electrolyte -

electrode head S7
cable length -
connector -

cell constant (cm⁻¹) -
cell type -

indicating metal platinum ring

Storage

InLab storage solution (30111142)
or 3 mol/L KCl (51350072)