


METTLER TOLEDO Use Case

Pharma - Compound Library




Customer

- Products for the treatment of lipid disorders and cardiovascular diseases
- Compound library with ~300'000 different substances



Application

- Preparing aliquots of 2-5mg, for (mostly) internal distribution
- 160-200 weighings per day



Need


- Fast weighing processes, stable results in no time
- Reliable and accurate results at minimal weights
- High user safety (substances potentially toxic)

Solution

- XP56 with low MinWeigh and high capacity
- Weighing directly into the tare containers
- ErgoClips micro for secure positioning of the tare containers and weighing directly, without any sample transfer
- MinWeigh Door micro for fast stabilization of the balance

Conclusion

- Save precious substances, due to one-step weighing process
- Save money (substance & time)
- Enhance productivity



Mettler-Toledo GmbH
Im Langacher
PO Box LabTec
8606 Greifensee
Switzerland

Subject to technical changes
©05/2006 Mettler-Toledo

www.mt.com/micro

For more information